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## Image Diagnosis Support Using Artificial Intelligence

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The AI (Artificial Intelligence) technologies are in progress for diagnostic support in radiology. A large quantity of radiological images continues being acquired by CT or MRI and so on. On the other hand, the radiologists, who read diagnostic images, are in short. And burdens to read images increase. And it is expected that diagnostic quality difference increases regionally. From this, expectation to diagnostic support by using AI including inspection efficiency and the reproducible improvement is high. Hitachi pushes forward development to make the whole inspection workflow of the radiology efficiency and high quality based on a concept called the DI × AI (Diagnostic Imaging with Artificial Intelligence) and is building the solution that utilized digital technology such as AI, Big data and IoT (Internet of Things). In this report, imaging position setting support of MRI and the reading diagnostic support of the lung cancer CT and brain MRI are shown.

**Key Word:** AI, Radiology, Diagnostic Support, Cancer, Brain

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- 1) Yokosawa S, et al.: Automated Scan Plane Planning for Brain MRI using 2D Scout Images. Proceedings of ISMRM, 3136, 2010.
- 2) Yokosawa S, et al.: Combining Rule-based and Machine Learning Approaches for Automated Scan Plane Planning of Spine MRI. JSRT, CyPos 10347, 2018.
- 3) 伊藤賢司、ほか: 拡散尖度画像と定量的磁化率画像を用いたパーキンソン症候群の早期鑑別診断. 第42回日本磁気共鳴医学会大会 (京都)、P-1-068, 2014.
- 4) Cabronero, J. A, et al.: In Vivo Quantitative Susceptibility Mapping (QSM) in Alzheimer's Disease. PLOS ONE, 8(11), e81093, 2013.
- 5) The National Lung Screening Trial Research Team: Reduced lung-cancer mortality with low-dose computed tomographic screening. N. Eng. J. Med., 365, 395-409, 2011.
- 6) U.S. Preventive Services Task Force home page, <https://www.uspreventiveservicestaskforce.org/Page/Document/UpdateSummaryFinal/lung-cancer-screening>
- 7) American College of Radiology home page, <https://www.acr.org/Clinical-Resources/Reporting-and-Data-Systems/Lung-Rads>
- 8) Nawa, T. et al.: A decrease in lung cancer mortality following the introduction of low-dose chest CT screening in Hitachi, Japan. Lung Cancer, 78 (3), 225-228, 2012.
- 9) Kusano, S. et al.: Efficacy of computer-aided diagnosis in lung cancer screening with low-dose spiral computed tomography: receiver operating characteristic analysis of radiologists' performance. Jpn. J. Radiol, 28 (9), 649-655, 2011.